# Sustainable compendial grade GMP detergent substitutes for Triton<sup>™</sup> X-100 in bioprocessing applications


Katherine McQuillan, Kritika Bajaj, Garima Sinha, Maria Broilo, Brian Malys, Miquel Mir, Sreejit Menon, Stephen Rumbelow, James Humphrey and Bradley Haltli **Croda International PIc & Associates** 

### Abstract

Triters "X 100 (coth) phone althoughts) has been used for decades as a detergent for virial isativisation in bioprocessing populations as the hexatine of human or niminal derived plasma, along with manuaccell lysis applications at laboratory scale and in cell-based biopharmecentrical production. While the ability of Tritors "X-100 (porform herea applications has never come indocusion, the environmental impact of this material automative levels in the ability of Tritors" X-100 (porform herea applications has never come indocusion, the environmental impact of this material automately led to its bain in Europe as of January 2021, with the exception of existing pharmaceutical applications allowed under a public exception or for small scale laboratory research user. This bain has less the significant ed under a public exemption or for small scale biodratory research user. This ban has led to significant in finding alternise detergents which are biodegradable. (MMP compliant and pharmaceutically plable for use in the manufacturing of new cell-derived drug products in Europe. Alternatives to Trion<sup>10</sup> Must perform equival in regard to call lysis and protein compatibility in order to ensure global applicability order to ensure global applicability wide range of bioprocessing applications. Employing xenotropic murmle leukemia virus (MruLV) and the catus PG4 cell line as a model lipic enveloped virus – host ystems, we compared the viral inactivation articles of earls PG4 cell line as a model lipic enveloped virus – host ystems, we compared the viral inactivation articles of earls High enveloped and early solid cell lines (PG4, PG4, PG4,


#### Introduction

Triton<sup>™</sup> X-100 is a 4-(1,1,3,3-tetramethylbuh)/johenyl ethoxylate detergent that is used widely in biopharmaceutical industry for virus inactivation and cell lysis applications. Triton<sup>™</sup> X-100 degrades to endoc disrupting by-products making if an aquatic reproductive toxin<sup>®</sup>. This led to its inclusion in the European Chem Reprovise candidation. Evaluation. Autorization of the Registrations. Evaluation. Autorization in Restriction of Chemicals (REACH) regulation<sup>®</sup>. Triton<sup>™</sup> X-100 was barned for use in the EU on Jan 4, 2021 a lew complions, requiring biopharmosatical annualizationes to final almostine detergents. Replacements Triton<sup>™</sup> X-100 have been developed<sup>®</sup> but there is room in the market to explore detrigents from "Crodia bord" triton<sup>™</sup> X-100 have been developed<sup>®</sup> but there is room in the market to explore detrigents from "Crodia bord" toxin and the centre of the centre of the regulation of the set of the centre of the regulation of the set of the centre of the tendent of the centre of th



or superior performance to Triton™ X-100, 31 o -100, 31 detergents belo ompared to Triton™ X-100 rgents, assays utilizing XI 25°. This study was cond their virus inactivation the virus inactivation used as a model for n properties and activity of the de r human retroviru To pr In the first phase we screened 31 detergents to identify the see we compared two prioritized detergents to industry ben ex TXR-1 and TXR-2 as exciting new products for virus inaction and biohomenet. performing detergen The results of this

## Viral inactivation testing



Virus inactivation was XMuLV. Tissue culture inf (TCID) was measured to viable viral particles after f the detergent. The reduction titer achieved by detergent with vira was achieved by detergent treatment was ated by determining the log reduction (LRF). A LRF of 3 was used for the screening assay. A LRF of 4 is an table target for detergent based inactivation. However, LRF of 8 - 10


0% 90% 99.9% 99.99% 99.999%

99.9999% 99.99999% 99.999999%

00 0000000%


## Detergent screening and potency testing


etergent materials belonging to 11 different chemical classes were screened for their viral inactivat Two chemically distinct products (Virodex™ TXR-1 and Virodex™ TXR-2) that exhibited the sa tivity as Triton™ X-100 were identif me concentration (0.025%) as Tritor 1. TXR-1 and TXR-2 achieved maximum vi X-100 after a 2 h detergent treatment at 22



#### Kinetics of viral inactivation

The time required to achieve a LRF of 4 is an important characteristic of a virus inactivation biomanufacturing process. A LRF of 4 is the industry standard virus inactivation target. Virus inactivation target. Virus inactivation process. A LRF of 4 is the industry standard virus inactivation process. A large of 4 is the industry standard virus inactivation process. A large of the virus inactivation process. A large of 4 is the industry standard virus inactivation industry and virus inactivation there is the virus inactivation to the virus inactivation that Triton <sup>TM</sup> - 100 after a 15 min exposure, with Virodex<sup>TM</sup> deterger an LRF greater than 3 after 3 5-minute treatment time, and after 60 min, the LRF increased to 6-8, exinitary standard target of 4. .... enect of temperature on showed equivalent or better Jex™ TXR-1 exhibiting more oth Viroder






n<sup>™</sup> X-100 is used to lyse cells uction cell lines. For this reason letter of una vecuba tutti productual tell inter- nu tans tessori, tee effectiveness or whode A efficiency of the second cell ly .... TOF CE Iuman Emi of con..... in the cell r

Detergent	CH	IO-K1	HEK-293T	
Virodex™ TXR-1	0.	0156	0.0156	
Virodex™ TXR-2	0.	0156	0.0625	
Triton™ X-100 control	0.	0156	0.0625	
CHO-K1 Cell Lysis		HEK-293T Cell Lysis		
Virodex <sup>16</sup> TXR-1 0.0039%	0.0156%	Virodex <sup>™</sup> TXR-1 0.03395	0.0156%	
Virodex™ TXR-2 0.0039%	0.0156%	Virodex™ TXR-2 0.0156%	0.0625%	


# Triton™ X-100 Triton™ X-100


# Protein compatibility of Virodex<sup>™</sup> TXR-1 and TXR-2

Biopharmaceutical bioprocesses are frequently used to manufacture protein-based products such as enzymes, protein therapeutics, growth factors and immunoglobins. Detergent treatments used to lyoe cells or nacivate lipid-devoloped viruses must not denature protein compatibility of Virodex<sup>10</sup> TKR-1 and TKR-2 was easing their effect on the activity of the screted enthypoinci alkaline phosphalase (SEAP) enzyme using a HEX-293 reporter cell line. Virodex<sup>10</sup> TKR-1 and TKR-2 was equivalent to Trindo-X-100. This is consistent with the non-toin nature of Virodex<sup>10</sup> determine deviated by the screte enthypoinci advance screter with the non-toin nature of Virodex<sup>10</sup> determine comparability equivalent to Trindo-X-100. This is consistent with the non-toin nature of Virodex<sup>10</sup> determine comparability equivalent to Trindo-X-100. This is consistent with the non-toin nature of Virodex<sup>10</sup> determine comparability explained to the color of 20% of the screter enthypoint cells and the protein comparability of 20% of the screte enthypoint cells and the protein comparability of 20% of the screte enthypoint cells and the protein comparability of 20% of the screte enthypoint cells and the constraints on point cells and the protein comparability of the cells and the cell of the screte enthypoint cells and the cell of the cells of the screte enthypoint cells and the cells of the (SEA), at concentra with the non-\*\* - SDS) sho




## Virodex<sup>™</sup> detergents outperform competitor products






ettors 2 and 3 were tested in separate head-to-head kinetic virus inactivation experiments with ex<sup>™</sup> detergents and Trion<sup>™</sup> X-100. Due to potent cytotoxicity against the *F*. catus PG4 cell line usec rus inactivation assay. Competitor 2 could not be tested at 0.1% and was instead tested at the non-to-ritation of 0.1%, in both comparisons, the **Virodex**<sup>™</sup> detergents out-performed the competitors.

## Performance consistency




Virodex™ removal: Protein A column affinity

In the manufacture of biotherapeutic proteins, it is typical to use chromatographic methods to purify the drug substance. A common reain used in the downstream processing of mAbs is the Protein A resin, with antboOs specificity. In mAb antrafustrue utiliary finition<sup>11</sup> × 10.01, tild detergent is waterial from the drug product alongaids other water materials prior to mAb elution. **Virodex<sup>111</sup> TXR-1** and **Virodex<sup>111</sup> TXR-2** were tested for their affinities. Nether **Virodex<sup>111</sup> TXR-1** or **Virodex<sup>111</sup> TXR-2** were found to have any affinity to the Protein A resin, and thus, the resin carb be effectively used to remove **Virodex<sup>111</sup> TXR-2** were found to have any affinity to the Protein A resin, and thus, the resin carb be effectively used to remove **Virodex<sup>111</sup> TXR-2**.




## Virodex<sup>™</sup> analytical detection

d sensitive detection and quantification methods are essential to confirm detergent removal cestical finished products. Croda developed a chromatography method that detects and TXR-1 and TXR-2 al low parts per taillion (pdp) concentrations using standard reversed-pl ure chromatography and mass spectrometry (UHPLC-MS) methodology. dex<sup>™</sup> TXR-1 the method achieved a 5 ppb limit of detection (LOD), a 10 ppb limit of quantification (LOQ) ear quantification response between 5 ppb and 100 ppb (R2 - 0.9998).



The same method achieved a 2 ppb LOD, a 5 ppb LOQ and a li 2 ppb – 1000 ppb (R<sup>2</sup> - 0.9992) for Virodex™ TXR-2.



## Freeze-thaw stability of Virodex™

Slability under freeze-thave conditions is essential for hopprocessing detergents, which may be subjected to repeated cycles of freezing and thaving during storage, transportation, or processing. The below data shows that **Wrodex** we detergents remain the execut-have stable after 14 cycles of freezing and thaving, and thus preserve their interprity under these conditions. During each cycle, solutions were frozen at a temperature of -20°C for 12 hours and thaved at 40°C for 12 hours.

Solutions tested	Results after 14 cycles	
Virodex™ TXR-1, 10.0% w/w	Clear solution, unchanged after 14 cycles	
Virodex™ TXR-1, 1.0% w/w	Clear solution, unchanged after 14 cycles	
Virodex™ TXR-2, 10.0% w/w	Clear solution, unchanged after 14 cycles	
Virodex™ TXR-2, 1.0% w/w	Hazy solution, unchanged after 14 cycles	

## Conclusion

n 16 ti

Of the 31 detergents screened for viral inactivation, two detergents (Virodex<sup>™</sup> TXR-1 and TXR-2) exhibited equivalent or superior performance to Triton<sup>™</sup> X-100 in all head-to-head comparisons. This data is consistent with previous research evaluating Triton™ X-100 alternatives<sup>6</sup>. Virodex™ detergents also outperformed competitor products in terms of safety and virus inactivation. Virodex™ detergents are effectively removed in mAb of satety and virus inactivation. Virodex<sup>11</sup> detergents are effectively removed in mAD manufacture using traditional Protein A chromatography. Sensitive ppb-level detection and quantification methods have been developed for industry standard instrumentation to aid customers in quality control analyses of finished products. Virodex<sup>10</sup> detergents also show no sign of degradation after repeated freeze-thaw cycles. The data presented here illustrates that the Virodex<sup>10</sup> detergents are exceptional alternatives to Triton<sup>10</sup> X-100 for backetmenoutling and hierarecencing and inductions. biopha rmaceutical and bioprocessing applications.

## References

122) ECHA. Available at: https://www.echa.europa.eu Available at: https://echa.europa.eu/substance-inform

- Hunter, A. K. et al. (20 Prog. 38: e3235 Created with BioRender.com



